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  
SAY SKIN

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A new substance,  
to be a mega-trend in future,

**PLASMA**

- Promotion of skin regeneration
- Strong sterilization effect
- Promotion of cosmetic absorption

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A mega-trend of beautifying  
apparatus,

**MICROCURRENT**

- Skin elasticity and pore reduction
- Promotion of physical metabolism

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A safe and eco-friendly  
ceramic,


**ZIRCONIA**

- A high-strengthened, precise and eco-friendly ceramic to emphasize safety

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

*Source : Non-Thermal Plasma for Acne and Aesthetic Skin Improvement*

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

- Clinical Institution : H&Bio SRC
- Clinical Period : June 7 ~ 21, 2021 (2 weeks)
- Test subject : Dermalab Full Skin 3cm x 3cm / Acety-Hexapeptide-8(Lys-FITC)
- Evaluation Item : The ability of absorption cosmetics active ingredients

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- Clinical Institution : H&Bio SRC
- Clinical Period : 2 weeks
- Test subject : Dermalab Full Skin 3cm x 3cm / Acety-Hexapeptide-8(Lys-FITC)
- Evaluation Item : The ability of absorption cosmetics active ingredients

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








- Clinical institution : OATC
- Clinical period : 2 weeks
- Clinical subject : 22 people, women aged 30~60 years old
- Evaluation item : The effect of skin elasticity for AURORA-lite with the ampoule

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
- Clinical institution : Korea Institute of Dermatological Sciences
- Clinical period : 6 weeks
- Test subject : 21 people, women aged 40~69 years old
- Evaluation item : Lifting of eye bag area under eyes, nasolabial area lifting, eye lifting, skin elasticity, skin density, skin tone

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
Before use    Immediately after use    6 weeks after use

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

Before use    Immediately after use    6 weeks after use

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Before use    Immediately after use    6 weeks after use

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




Please click the image below.

1. After cleansing your face, remove moisture on your face with a towel.

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© 2021 WITHNIX. All rights reserved
All intellectual property rights, including copyrights to the above document, are at Withnix Co., Ltd.