


SAY SKIN

AURORA-me

Quick Guide to using the device AURORA-me

HOW TO USE

[Switch it On/Off]



Press On/Off & Level button for 2 seconds

[Choose Level]

Press On/Off & Level button shortly after the power has been turned on.

Each time you press the button, the level will go up. If you press the button at Level 3, it will go back to Level 1.

- Level 1 (White LED)
- Level 2 (Orange LED)
- Level 3 (Red LED)

STEP 01.



Tone your facial skin after cleansing.

STEP 02.



Apply SAY SKIN super moisture cream on your skin.
(Or you may use other skin care products. Applying cream will reduce friction on your skin while using the device and promote the absorption of effective ingredients of skin care products.)

STEP 03.



Gently put the device on the area you want and move the device as if you scanned your face with it following your skin texture from bottom to top.